



Quinolones Rapid Test Kit

(Art. NO. STW006.18.01)

【Product description】

This Kit (Art. NO. STW006.18.01) is used for the qualitative determination of Quinolones in **muscle tissue** based on the colloidal gold immunochromatography technology. The test has been validated for Chicken,Pork,Fish,Shrimp,Duck . The Quinolones rapid test kit is convenient, fast, sensitive, and can be used for on-site large-scale screening.

【Detection Limit】

Chicken,Pork,Fish,Shrimp,Duck:

Quinolones	Detection limit
Enrofloxacin	40 ppb (ng/g)
Ciprofloxacin	40 ppb (ng/g)
Norfloxacin	40 ppb (ng/g)
Flumequine	40 ppb (ng/g)
Lomefloxacin	40 ppb (ng/g)
Danofloxacin	40 ppb (ng/g)
Oxolinic acid	40 ppb (ng/g)
Enoxacin	40 ppb (ng/g)
Pefloxacin	40 ppb (ng/g)
Marbofloxacin	40 ppb (ng/g)
Ofloxacin	40 ppb (ng/g)

【Components provided】

1. Test strip, 40 pcs / kit.
2. Sample extractant A
3. Sample extractant B
4. reconstitution fluid
5. Instruction manual

【Storage instructions】

1. Store the kit at 20±5°C (59~77°F), and do not freeze the components of the test kit.



2. The test is valid for 12 months. No quality guarantee can be assumed after the expiration date on the kit label.

【Materials Required but not Provided】

1. Equipment:

- 1.1 Graduated micropipette / tips 20-200 μ L
- 1.2 Electronic balance (accuracy: 0.01 g)
- 1.3 Vortex mixer
- 1.4 Centrifuge (4000 g)
- 1.5 Nitrogen Evaporator
- 1.6 Timer

【Test Procedure】

1 Preparation


- 1.1 Please read the entire manual carefully before starting your experiment. Take out the test strips needed and the samples to reach room temperature ($20\pm 5^{\circ}\text{C}/59\sim 77^{\circ}\text{F}$).
- 1.2 Weigh 1 ± 0.05 g of homogenized sample into 10 mL centrifuge tube.
- 1.3 Add 0.8 mL sample extractant A, mix vigorously by vortex mixer for 1 min.
- 1.4 Add 4 mL sample extractant B, mix vigorously by vortex mixer for 3 min.
- 1.5 Then centrifuge the sample for 5 min with $4000\times g$ at room temperature.
- 1.6 Transfer 0.5 mL sample supernatant into a new 5 mL centrifuge tube and evaporate it with nitrogen stream in 60°C (140°F) water bath.
- 1.7 Add 0.8 mL reconstitution fluid, and then mix adequately by vortex mixer for 30 s.

2 Testing

- 2.1 Take out the quantity needed of strips and label them, add 120 μ L sample into the sample well (S).
- 2.2 5 min after sample addition, interpret the result according to the color intensities of C Line (Control Line) and T Line (Test Line). (See Result Determination).

【Result interpretation】

1. Visual interpretation



Negative: T Line and C Line are both red, and color of T Line is deeper than C Line or the same with C Line.

Positive: C Line is red, T line does not show any color or the color is obviously lighter than C Line.

Invalid: Strip does not develop a C Line. The corresponding sample should be re-tested.

2. Interpretation with strip reader

Please refer to the manual of reader.

【Notes】

1. No water or other liquid can be used as negative control.
2. The sample should be homogenized to facilitate the extraction of the drug.
3. Test should be finished within 1 hour.
4. This strip is for single use only.
5. Please do not use if the package is broken, or dirty or out of expiration date.
6. This kit is for screening test only, positive result should be further confirmed with other methods.
7. All wastes should be properly decontaminated prior to disposal. Disposal should be in accordance with local, regional, and national regulations.

WDWK makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, WDWK will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. WDWK shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

【Manufacturer】

Beijing WDWK Biotechnology Co., Ltd

Add: Building 3, Courtyard 9, Dijin Road, Haidian District, Beijing, China

Zip code: 100095

Tel: +86-10-62974201



Fax: +86-10-62987854

E-mail: info@wdwkbio.com

Website: www.wdwkbio.com